

TURF TOE

YOU WILL NEED





STRIP ONE: Tear backing paper two inches from end to create starting anchor.

STRIP TWO: Tear backing paper at the center of the half strip to begin application.

ALL COLORS of KT Tape® are designed to perform equally.



TO CREATE A HALF STRIP: Fold a full strip in half and cut folded edge rounding the corners.

BEFORE YOU START

APPLY BEFORE ACTIVITY Apply one hour before beginning activity

> **CLEAN SKIN** Clean dirt, oils and lotions from area with rubbing alcohol

ACTIVATE ADHESIVE

After application rub tape vigorously to activate adhesive





While sitting, cross your ankle over your opposite knee.

PREPARE:

Cut a strip of tape vertically to create two, one-inch wide strips.



ANCHOR: Apply anchor on the bottom of your big toe with no stretch.

STRIP TWO



ANCHOR: Tear backing paper of second strip in its center. Guide it between your first and second toes with a



stretch, apply up the

bottom of the foot and

heel with no stretch on

along the back of the

the opposite end.

APPLY: Next lay day the ends underneath the toe digonally with no stretch.





CAUTION: If you have skin sensitivities, cancer, or are pregnant, consult your doctor before use. Discontinue use if skin becomes irritated or sore. K TAPE® is not a replacement for professional medical care. Warranties and remedies limited to product replacement cost. READ ALL CAUTIONS ON ENCLOSED INSTRUCTION SHEET PRIOR TO USE.